



Acquire


napari workshop on multi-dimensional optical microscopy

Alan Liddell, Chan Zuckerberg Initiative

2023-09-27

CZI Imaging Tech

How CZI is involved in the napari Community



napari viewer
visualization & analysis interface

CZI: Supports some core
developers and contributes to repo

napari hub
discovery & sharing
of image analysis plugins


CZI: Develops in service to the
community and **funds plugin grants**



Give researchers
access to
reproducible,
quantitative image
analysis.

Acquire

- Exploratory project started in 2022
- 2-camera simultaneous streaming
- High speed camera streaming (2 GiB / s)
- Cloud-optimized large data storage (Zarr)




Collaboration with
[Allen Institute](#)
& [SF Biohub](#)

<https://github.com/acquire-project/acquire-python>



Acquire + napari

Leverage the viewer and widget support



Live streaming
to napari

Instrument
parameters


Laser control

Sample stage
location





Acquire + napari

Simultaneous two-camera acquisition



Acquire architecture



Using Acquire in Python, part 1



```
import acquire

runtime = acquire.Runtime()
dm = runtime.device_manager()
config = runtime.get_configuration()

# configure stream 0
stream = 0

config.video[stream].camera.identifier = dm.select(acquire.DeviceKind.Camera, "simulated: radial
sin")
config.video[stream].camera.settings.binning = 1
config.video[stream].camera.settings.exposure_time_us = 5e4
config.video[stream].camera.settings.pixel_type = acquire.SampleType.U8
config.video[stream].camera.settings.shape = (1024, 768)
config.video[stream].storage.identifier = dm.select(acquire.DeviceKind.Storage, "Zarr")
config.video[stream].storage.settings.filename = "output.zarr"
config.video[stream].storage.settings.chunking.max_bytes_per_chunk = 32 * 2**20
config.video[stream].max_frame_count = 100
```

Using Acquire in Python, part 2



```
import logging

# start the runtime
runtime.start()
print(runtime.get_state()) # DeviceState.Running

# inspect and work with the data
if a := runtime.get_available_data(stream):
    packet = a.get_frame_count()
    for f in a.frames():
        frame_data = f.data() # frame_data is a numpy array
        # do work with frame_data
        del f # manually release frame
    del a # manually release available data

runtime.stop()
```


Acquire with μ Manager


Enable a wider user base to access Acquire



Micromanager

Widely used open-source, Java-based image acquisition software

1. Simultaneous 2 camera acquisition
2. Access to existing device control
3. Enables streaming to OME-Zarr within μ Manager



Acknowledgement

CZI Imaging tech team
Allen Institute for Neural Dynamics
Politecnico Milano

Learn more about Acquire

[Acquire Github repos](#)

ExA-SPIM in [Nature](#) ([preprint](#))